



Proactive Health

SINGLE-GENE RISK™

84 GENES

TEST CODE: PR22020

Overview

MyOme Proactive Health Single-Gene risk report uses a Blended Genome-Exome (BGE) backbone built from whole exome sequencing and low coverage whole genome sequencing to identify a range of variant types. This allows MyOme to re-query a patient's data as healthcare needs change and new information about the genome is discovered.

Clinical Use

Test is intended for a wellness screening of germline heritable conditions in individuals from an asymptomatic population. MyOme annotates and interprets variants according to American College of Medical Genetics (ACMG) guidelines¹, and reports pathogenic or likely pathogenic variants. Genetic testing may provide information about an individual's disease risk, support a clinical diagnosis, and assist with the development of a personalized treatment and management strategy in conjunction with standard clinical assessment.

Method

PCR-free whole genome library is constructed and a sub-aliquot is taken through PCR amplification and exome selection. The blended genome and exome libraries are sequenced to generate 2x150 bp paired-end reads resulting in low-coverage whole genome and higher coverage exome data. In-house pipeline allows identification of single-nucleotide variants (SNVs) and small insertions and deletions (indels). Variant interpretation by qualified scientists based on guidelines by the ACMG.

Sample Types

- Blood (2 EDTA tubes)
- Saliva (2 tubes)
- Buccal (3 swabs)

Turn Around Time

- From initial sample received, approximately 6 to 8 weeks
- For previously processed sample, approximately 2 to 4 weeks

Included

- Analysis of SNVs and small insertions and deletions
- Confirmation of Pathogenic/Likely Pathogenic variants by orthogonal method (e.g., Sanger sequencing)
- Cohesive report with actionable recommendations
- **84 Genes included:** *ACTA2, ACTC1, ACVRL1, APC, APOB, ATM, ATP7B, BAG3, BMPR1A, BRCA1, BRCA2, BTBD, CACNA1S, CALM1, CALM2, CALM3, CASQ2, CHEK2, COL3A1, DES, DSC2, DSG2, DSP, ENG, FBN1, FLNC, GAA, GLA, HFE, HNF1A, KCNH2, KCNQ1, LDLR, LDLRAP1, LMNA, MAX, MEN1, MLH1, MSH2, MSH6, MUTYH, MYBPC3, MYH11, MYH7, MYL2, MYL3, NF2, OTC, PALB2, PCSK9, PKP2, PMS2, PRKAG2, PTEN, RB1, RBM20, RET, RPE65, RYR1, RYR2, SCN5A, SDHAF2, SDHB, SDHC, SDHD, SMAD3, SMAD4, STK11, TGFBF1, TGFBF2, TMEM127, TMEM43, TNNT3, TNNT1, TNNT2, TP53, TPM1, TRDN, TSC1, TSC2, TTN, TTR, VHL, WT1*

Test Performance²

- ≥60x average exome-wide coverage
- ≥1x average genome-wide coverage
- ≥90% of exonic regions at ≥20x depth
- >99.5% sensitivity for SNVs
- >97.5% sensitivity for small insertions/deletions

1. American College of Medical Genetics and Genomics. SF v3.2 list for reporting of secondary findings in clinical exome and genome sequencing: A policy statement of the ACMG. *Genet Med*. June 22, 2023. doi: 10.1016/j.gim.2023.100866.

2. MyOme Inc, Data on File.

The test described above was developed, and its performance characteristics determined, by MyOme, Inc., a clinical laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) to perform high complexity clinical laboratory testing. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA).

MyOme is not responsible for the content or accuracy of third-party websites.

CONDITION-GENE RELATIONSHIP

The genes listed below are analyzed in this report. MyOme selected them based on our Gene Inclusion Framework Guidelines. Genes are prioritized based on clinical validity, actionability, penetrance/prevalence, and feasibility.

Cardiovascular	
Condition	Gene(s)
Arrhythmogenic right ventricular cardiomyopathy	DES, DSC2, DSG2, DSP, PKP2, TMEM43
Brugada syndrome	SCN5A
Catecholaminergic polymorphic ventricular tachycardia	CALM1, CALM2, CALM3, CASQ2, RYR2, TRDN
Dilated cardiomyopathy	ACTC1, BAG3, DES, FLNC, LMNA, MYH7, SCN5A, TNNC1, TNNI3, TNNT2, TPM1, TTN, RBM20
Ehlers-Danlos syndrome, vascular type	COL3A1
Emery-Dreifuss muscular dystrophy	LMNA
Fabry disease	GLA
Familial hypercholesterolemia	APOB, LDLR, LDLRAP1, PCSK9
Familial thoracic aortic aneurysm and dissection	ACTA2, MYH11, SMAD3
Hereditary transthyretin-related amyloidosis	TTR
Hypertrophic cardiomyopathy	ACTC1, MYBPC3, MYH7, MYL2, MYL3, PRKAG2, TNNC1, TNNI3, TNNT2, TPM1
Loeys-Dietz syndrome	TGFBR1, TGFBR2, SMAD3
Long QT syndrome	CALM1, CALM2, CALM3, KCNH2, KCNQ1, SCN5A, TRDN
Marfan syndrome	FBN1
Myofibrillar myopathy	BAG3, DES, FLNC
Short QT syndrome	KCNH2, KCNQ1

Other	
Condition	Gene(s)
Biotinidase deficiency	BTD
Hereditary hemochromatosis	HFE
Hereditary hemorrhagic telangiectasia	ACVRL1, ENG, SMAD4
Malignant hyperthermia	CACNA1S, RYR1
Monogenic diabetes	HNF1A
Ornithine transcarbamylase deficiency	OTC
RPE65-related retinopathy	RPE65
Wilson disease	ATP7B

Cancer	
Condition	Gene(s)
Familial adenomatous polyposis	APC
Familial ovarian cancer	PALB2
Gastrointestinal stromal tumor	KIT
Hereditary breast cancer	ATM, CHEK2, PALB2
Hereditary breast and ovarian cancer	BRCA1, BRCA2
Hereditary nonpolyposis colon cancer	ATM
Hereditary paraganglioma-pheochromocytoma syndrome	MAX, SDHAF2, SDHB, SDHC, SDHD, TMEM127
Juvenile polyposis syndrome	BMPR1A
Juvenile polyposis with hereditary hemorrhagic telangiectasia	SMAD4
Li-Fraumeni syndrome	TP53
Lynch syndrome	MLH1, MSH2, MSH6, PMS2
Multiple endocrine neoplasia	MEN1, RET
MUTYH-associated polyposis	MUTYH
Neurofibromatosis type 2	NF2
Peutz-Jeghers syndrome	STK11
PTEN hamartoma tumor syndrome	PTEN
Retinoblastoma	RB1
Tuberous sclerosis complex	TSC1, TSC2
Von Hippel-Lindau syndrome	VHL
WT1-related Wilms tumor	WT1

The test described above was developed, and its performance characteristics determined, by MyOme, Inc., a clinical laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) to perform high complexity clinical laboratory testing. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA).
MyOme is not responsible for the content or accuracy of third-party websites.



TEST CODE: PR32005

Overview

MyOme Proactive Health Medication Response report uses a Blended Genome-Exome (BGE) backbone to identify variants associated with medication response. This allows MyOme to re-query a patient's data as healthcare needs change and new information about the genome is discovered.

Clinical Use

Test is intended to facilitate the use of pharmacogenomic guidance in a general care setting. MyOme reports variants and star alleles in 15 pharmacogenes that are known to impact an individual's response to medication. The results of this test should be interpreted by a trained healthcare provider based on the full context of a patient's medical situation.

Method

Genomic DNA obtained from submitted samples is used to construct a PCR-free whole genome library and an exome library and sequenced using Illumina technology. Reads are aligned to the GRCh38 reference assembly. MyOme's Medication Response variant calling pipeline is used to call variants, star alleles, copy number variants, and assign pharmacogenomic diplotypes.

Sample Types

- Blood (2 EDTA tubes)
- Saliva (2 tubes)
- Buccal (2 swabs)

Turn Around Time

- From initial sample received, approximately 6-8 weeks

Included

- Analysis of key genes with known pharmacogenomic implications for more than 70 drugs, including antidepressants, statins and opioids
- 96% Tier 1 variant alleles recommended by the Association for Molecular Pathology (AMP) PGx Working Group for *CYP2C19*¹, *CYP2C9*², *CYP2D6*³ (excluding deletions/duplications), *TPMT* and *NUDT15*⁴
- All Tier 2 variant alleles recommended by the AMP PGx Working Group for *CYP2C19*¹, *CYP2C9*², *CYP2D6* (excluding hybridizations)³
- Independent healthcare provider review of analysis to make any potential drug recommendations based on gene-drug interactions outlined in the Clinical Pharmacogenetics Implementation Consortium (CPIC)⁵ guidelines and FDA Table of Pharmacogenomic Associations⁶
- 15 genes analyzed: ***CYP2B6***: *4; *6; *9; *18; *22 ***CYP2C9***: *2; *3; *4; *5; *6; *8; *11; *12; *13; *15; *16; *26; *28; *29; *30; *31; *42; *55 ***CYP2C19***: *2; *3; *4; *5; *6; *7; *8; *9; *10; *17; *35 ***CYP2D6***: *2; *3; *4; *6; *7; *8; *9; *10; *11; *12; *14; *15; *17; *21; *31; *40; *41; *42; *49; *56; *59; *100; *114; ***CYP3A4***: *22 ***CYP3A5***: *3; *6 ***CYP4F2***: *3 ***DPYD***: rs3918290; rs55886062; rs59086055; rs67376798; rs75017182+rs56038477; rs112766203; rs115232898; rs146356975; rs183385770 ***F5***: rs6025 ***IFNL3***: rs12979860 ***NUDT15***: *3; *4; *9, ***SLCO1B1***: *5; *9; *14; *20 ***TPMT***: *2; *3A; *3B; *3C; *4; *11; *29 ***UGT1A1***: *6; *27 ***VKORC1***: rs9923231

Test Performance⁷

- >98.5% accuracy for diplotypes
- 98% of core pharmacogenomic variants at ≥20x depth
- 60x average exome-wide coverage

1. Pratt VM et al. Recommendations for Clinical CYP2C19 Genotyping Allele Selection: A Report of the Association for Molecular Pathology. J Mol Diagn. 2018 May;20(3):269-276. doi: 10.1016/j.jmoldx.2018.01.011. 2. Pratt VM et al. Recommendations for Clinical CYP2C9 Genotyping Allele Selection. J Mol Diagn. 2019 Sep;21(5):746-755. doi: 10.1016/j.jmoldx.2019.04.003. 3. Pratt VM et al. Recommendations for Clinical CYP2D6 Genotyping Allele Selection. J Mol Diagn. 2021 Sep;23(9):1047-1064. doi: 10.1016/j.jmoldx.2021.05.013. 4. Pratt VM et al. TPMT and NUDT15 Genotyping Recommendations. J Mol Diagn. 2022 Oct;24(10):1051-1063. doi: 10.1016/j.jmoldx.2022.06.007. 5. CPIC. Genes-Drugs. Web. cpicpgx.org/genes-drugs. Accessed 2024 Dec. 6. US Food & Drug Administration. Table of Pharmacogenetic Associations. 2022 Oct. Web. Accessed 2024 Dec. 7. MyOme Inc, Data on File

The test described above was developed, and its performance characteristics determined, by MyOme, Inc., a clinical laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) to perform high complexity clinical laboratory testing. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA).

MyOme is not responsible for the content or accuracy of third-party websites.

DRUG-GENE IMPACT

The list of drugs below can lead to a major or moderate drug-gene interaction based on guidelines described above.^{5,6}

Behavioral Health	
Drug Name	Gene(s)
amitriptyline	CYP2D6, CYP2C19
amoxapine	CYP2D6
amphetamine	CYP2D6
aripiprazole	CYP2D6
aripiprazole lauroxil	CYP2D6
atomoxetine	CYP2D6
brexpiprazole	CYP2D6
citalopram	CYP2C19, CYP2D6
clobazam	CYP2C19
clomipramine	CYP2D6
clozapine	CYP2D6
desipramine	CYP2D6
diazepam	CYP2C19
doxepin	CYP2D6, CYP2C19
duloxetine	CYP2D6
escitalopram	CYP2C19
fluoxetine	CYP2D6
fluvoxamine	CYP2D6
iloperidone	CYP2D6
imipramine	CYP2D6
lofexidine	CYP2D6
nortriptyline	CYP2D6
paroxetine	CYP2D6
perphenazine	CYP2D6
protriptyline	CYP2D6
sertraline	CYP2C19
thioridazine	CYP2D6
trimipramine	CYP2D6, CYP2C19
valbenazine	CYP2D6
venlafaxine	CYP2D6
vortioxetine	CYP2D6
Urology	
Drug Name	Gene(s)
darifenacin	CYP2D6
fesoterodine	CYP2D6
mirabegron	CYP2D6
tamsulosin	CYP2D6
tolterodine	CYP2D6
Reproductive and Sexual Health	
Drug Name	Gene(s)
flibanserin	CYP2C19, CYP2C9, CYP2D6

Gastroenterology	
Drug Name	Gene(s)
dexlansoprazole	CYP2C19
dronabinol	CYP2C9
esomeprazole	CYP2C19
lansoprazole	CYP2C19
metoclopramide	CYP2D6
nateglinide	CYP2C9
omeprazole	CYP2C19
ondansetron	CYP2D6
pantoprazole	CYP2C19
rabeprazole	CYP2C19
Infectious Disease	
Drug Name	Gene(s)
atazanavir	UGT1A1
dolutegravir	UGT1A1
efavirenz	CYP2B6
peginterferon alfa-2a	IFNL3
peginterferon alfa-2b	IFNL3
quinine	CYP2D6
voriconazole	CYP2C19
Cardiology	
Drug Name	Gene(s)
atorvastatin	SLC01B1
carvedilol	CYP2D6
clopidogrel	CYP2C19
fluvastatin	CYP2C9, SLC01B1
lovastatin	SLC01B1
mavacamten	CYP2C19
pitavastatin	SLC01B1
pravastatin	SLC01B1
rosuvastatin	SLC01B1
simvastatin	SLC01B1
warfarin	CYP2C9, CYP4F2, VKORC1
Miscellaneous	
Drug Name	Gene(s)
abrocitinib	CYP2C19
cevimeline	CYP2D6
eliglustat	CYP2D6
lesinurad	CYP2C9
meclizine	CYP2D6
tropisetron	CYP2D6

Neurology	
Drug Name	Gene(s)
brivaracetam	CYP2C19
deutetrabenazine	CYP2D6
donepezil	CYP2D6
fosphenytoin	CYP2C9
modafinil	CYP2D6
phenytoin	CYP2C9
pimozide	CYP2D6
pitolisant	CYP2D6
siponimod	CYP2C9
tetrabenazine	CYP2D6
valbenazine	CYP2D6
Pain Management	
Drug Name	Gene(s)
carisoprodol	CYP2C19
celecoxib	CYP2C9
codeine	CYP2D6
elagolix	SLC01B1
flurbiprofen	CYP2C9
ibuprofen	CYP2C9
lornoxicam	CYP2C9
meloxicam	CYP2C9
oliceridine	CYP2D6
piroxicam	CYP2C9
tenoxicam	CYP2C9
tramadol	CYP2D6
Hematology/Oncology	
Drug Name	Gene(s)
capecitabine	DPYD
erdafitinib	CYP2C9
eltrombopag	F5
fluorouracil	DPYD
gefitinib	CYP2D6
mercaptopurine	TPMT, NUDT15
tamoxifen	CYP2D6
thioguanine	TPMT, NUDT15
Transplant	
Drug Name	Gene(s)
azathioprine	TPMT, NUDT15
tacrolimus	CYP3A5

5. CPIC. Genes-Drugs. Web. cpicpgx.org/genes-drugs. Accessed 2024 Dec. 6. US Food & Drug Administration. Table of Pharmacogenetic Associations. 2022 Oct. Web. Accessed 2024 Dec. The test described above was developed, and its performance characteristics determined, by MyOme, Inc., a clinical laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) to perform high complexity clinical laboratory testing. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA). MyOme is not responsible for the content or accuracy of third-party websites.

Jane Doe

Biological Sex: Female
Date of Birth: 11/13/2018
Sample ID: SM12805
Sample Type: BLOOD
Collection Date: 11/13/2018
Received Date: 11/14/2018

Clinic: The City Clinic

Physician: Jane Smith, M.D.
Phone: 555-555-5555
Fax: 555-555-5555
NPI: 0123456789

Requisition ID:

RQ12345

Report Number:

RP12345

Report Date:

11/01/2024

Scan to learn more



<https://myome.com/r/?id=RP12345>

Your lifetime risk of developing breast cancer is average (less than 20%) based on your integrated score.

CLINICAL CONTEXT: This test integrates known clinical risk factors and a polygenic risk score. It does NOT incorporate single gene findings in breast cancer predisposition genes.

RISK DETAILS

	LIFETIME RISK	5-YEAR RISK
Integrated Risk	14.2%	0.5%
Clinical Risk	12.2%	0.4%
General Population Risk	11.7%	0.1%

Integrated Risk: The risk of developing breast cancer based on the combination of a polygenic risk score and the Tyrer-Cuzick clinical risk model.

Clinical Risk: The risk of developing breast cancer based on the Tyrer-Cuzick clinical risk model.

General Population Risk: The average risk of developing breast cancer for a biological female in the general population of the same age.

NEXT STEPS

- Next steps described in this section are based on lifetime risk as estimated by the integrated risk score.
- Individuals with an estimated lifetime breast cancer risk less than 20% can typically follow the general population recommendations from society guidelines, as outlined below. Please discuss the risk score as part of a comprehensive risk assessment with a healthcare provider, as some factors such as prior chest wall radiation, or mammographic breast density, may warrant additional breast cancer screening [1,4,5].
- Option to start annual mammograms at age 40 with recommendation to start by age 45.
- These results should be interpreted in the context of the individual's personal medical history or family history.

REFERENCES

1. Saslow et al. *American Cancer Society guidelines for breast screening with MRI as an adjunct to mammography*. CA Cancer J Clin. 2007; 57(57):75-89. PMID: 17392385.
2. Tshiaba et al. *Integration of a Cross-Ancestry Polygenic Model With Clinical Risk Factors Improves Breast Cancer Risk Stratification*. JCO Precis Oncol. 2023; 7(7):e2200447. PMID: 36809055.
3. Tyrer et al. *A breast cancer prediction model incorporating familial and personal risk factors*. Stat Med. 2004; 23(23):1111-30. PMID: 15057881.
4. Oeffinger et al. *Breast Cancer Screening for Women at Average Risk: 2015 Guideline Update From the American Cancer Society*. JAMA. 2015; 314(314):1599-614. PMID: 26501536.
5. US et al. *Screening for Breast Cancer: US Preventive Services Task Force Recommendation Statement*. JAMA. 2024; 331(331):1918-1930. PMID: 38687503.
6. Monticciolo et al. *Breast Cancer Screening for Women at Higher-Than-Average Risk: Updated Recommendations From the ACR*. J Am Coll Radiol. 2023; 20(20):902-914. PMID: 37150275.

TEST METHODS

- Patient data is provided by the ordering physician. Specimen receipt, accessioning, data analysis and interpretation is performed by MyOme, Inc., 1455 Adams Dr., Ste 1150, Menlo Park, CA 94025, CLIA#05D2203070. Whole Genome Sequencing, excluding data analysis and interpretation, is performed by Broad Clinical Labs LLC, 27 Blue Sky Dr, Burlington, MA 01803, CLIA#22D2055652.
- Genomic DNA obtained from submitted samples was sequenced using Illumina technology. Reads were aligned to the NCBI GRCh37.p13 reference assembly.
- A polygenic risk score (PRS) is calculated for each of 5 continental ancestries of which the patient is a part as the sum of the patient's risk alleles weighted by the allele-specific effect sizes. The raw scores are centered using four principal components and standardized with a population-specific standard deviation. Standardized PRSs weighted by fractional ancestry and ancestry-specific effect sizes are summed (caPRS) [2].
- The standardized PRS is integrated with the risk based on the Tyrer-Cuzick (TC) model to estimate a 5-year and remaining lifetime risk of developing breast cancer [3].
- **This tool cannot be used to detect rare pathogenic variants including those in hereditary cancer predisposition genes.**

TEST LIMITATIONS

- The results of this test may not be valid if the patient has a pathogenic variant in a breast cancer predisposition gene. The integrated risk score estimate does not account for pathogenic variants in genes with limited or disputed breast cancer association (e.g. BRIP1), as the available data is currently insufficient to accurately quantify the breast cancer risk associated with these genes.
- The clinical risk based on the Tyrer-Cuzick risk model was calculated based on the patient data provided by the ordering physician. Incorrect or missing information will impact this calculation and the integrated risk score.
- A risk calculation will not be performed for biological males. A risk calculation will also not be performed for biological females who are under the age of 18, over the age of 84, known to carry a pathogenic variant in a breast cancer predisposition gene or have a personal history of breast cancer.
- A risk calculation will not be performed when there is missing information necessary to perform the calculation, including but not limited to age and first degree family history of breast cancer.
- The breast cancer integrated risk is a risk assessment tool NOT a diagnostic. These results should be interpreted in the context of the individual's personal medical history and family history.
- Performance of this tool may be reduced in certain populations.
- Like most tests, this test carries a risk of false negative or false positive results, which may be caused by, without limitation, sample contamination from biological or non-biological sources, specimen marking issues, rare genetic variants interfering with analysis, and other technical issues and limitations.

DISCLAIMERS

This test was developed, and its performance characteristics were determined, by MyOme, Inc., a clinical laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) to perform high complexity clinical laboratory testing. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA).

George Washington Carver

Sex Assigned at Birth: Male
Date of Birth: 06/15/1967
Sample ID: SM-0000232
Sample Type: BLOOD
Collection Date: 02/01/2024
Received Date: 02/03/2024

Clinic: Cardiovascular Health Center

Physician: Maria Chen, M.D., Ph.D.
Phone: 510-555-0000
NPI: 1234567890

Requisition ID:

RQ0000059

Report Number:

RP12345

Report Date:

02/13/2025

Scan to learn more




<https://myome.com/r/?id=RP12345>

High Risk

Based on the integrated risk score, this patient has a 31.8% chance of experiencing a coronary artery disease (CAD) event in the next 10 years.

10-YEAR ABSOLUTE RISK OF CAD



Integrated Risk: The probability of having a CAD-related event within the next 10 years based on the combination of a polygenic risk score (PRS) and clinical risk factors.

Clinical Risk: The probability of having a first atherosclerotic cardiovascular disease (ASCVD) event within the next 10 years based on clinical risk factors. Factors include age, race, total and HDL cholesterol, systolic blood pressure, smoking and diabetes status, and whether the patient is receiving treatment for blood pressure abnormalities.

Based on the integrated risk score, this patient has a 31.8% chance of experiencing a coronary artery disease (CAD) event in the next 10 years.

This test determined that the patient is at high risk for CAD, defined as an absolute 10-year integrated risk greater than 20%. Being at high risk for CAD does not mean that the patient will definitely develop the disease. There may be other genetic and non-genetic factors not considered here that influence the patient's risk.

NEXT STEPS

The AHA recommends all patients follow a heart-healthy lifestyle, including smoking cessation, regular physical activity, and a balanced diet [3].

For patients with a high 10-year ASCVD risk estimation, the AHA recommends starting a high-intensity statin and addressing comorbid conditions such as hypertension or diabetes [3].

CLINICAL MEASUREMENTS

The clinical measurements used to estimate the patient's risk of developing CAD were provided by the ordering physician and are detailed below. Changes in these measurements will result in a change in the patient's risk estimation.

Age	57	Total Cholesterol	131 mg/dL
Sex Assigned at Birth	Male	HDL Cholesterol	39 mg/dL
Self-reported Race	Black/African American	History of Diabetes	Yes
Systolic Blood Pressure	118 mm Hg	Smoking Status	Current smoker
Diastolic Blood Pressure	76 mm Hg	Hypertension Treatment	No

REFERENCES

1. Tshiaba, Placede et al. Genetics in Medicine Open, Volume 1, Issue 1, 100361
2. Goff et al. *2013 ACC/AHA guideline on the assessment of cardiovascular risk: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines*. Circulation. 2014; 129(129):S49-73. PMID: 24222018.
3. Arnett et al. *2019 ACC/AHA Guideline on the Primary Prevention of Cardiovascular Disease: Executive Summary: A Report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines*. Circulation. 2019; 140(140):e563-e595. PMID: 30879339.

TEST METHODS

- Patient data is provided by the ordering physician. Specimen receipt, accessioning, data analysis and interpretation is performed by MyOme, Inc., 1455 Adams Dr., Ste 1150, Menlo Park, CA 94025, CLIA#05D2203070. Whole Genome Sequencing, excluding data analysis and interpretation, is performed by Broad Clinical Labs LLC, 27 Blue Sky Dr, Burlington, MA 01803, CLIA#22D2055652.
- Genomic DNA obtained from submitted samples is sequenced using Illumina technology. Reads are aligned to the NCBI GRCh38 reference assembly.
- A polygenic risk score (PRS) is calculated for each of 5 continental ancestries of which the patient is a part as the sum of the patient's risk alleles weighted by the allele-specific effect sizes. The raw scores are centered using four principal components and standardized with a population-specific standard deviation. Standardized PRSs weighted by fractional ancestry and ancestry-specific effect sizes are summed (caPRS) [1].
- The standardized caPRS is integrated with the clinical risk based on the Atherosclerotic cardiovascular disease (ASCVD) Pooled Cohort Equations (PCE) model [2] to estimate the patient's 10-year risk of developing CAD.

TEST LIMITATIONS

- The clinical risk based on the ASCVD PCE risk model was calculated using the patient data provided by the ordering physician. Incorrect information will impact this calculation and the integrated risk score.
- A risk calculation will not be performed for individuals younger than 40 or older than 79 years old or for individuals missing clinical measurements or individuals with a personal history of CAD.
- This test is a risk assessment tool NOT a diagnostic. These results should be interpreted in the context of the individual's personal medical history and family history.
- Performance of this tool may be reduced in certain populations.
- A history of stem cell or bone marrow transplantation, or recent blood transfusion may impact the accuracy of the results.
- Testing is unavailable for samples damaged by human error or lost/destroyed due to weather, transit issues or other problems beyond the control of MyOme.
- Like most tests, this test carries a risk of false negative or false positive results, which may be caused by, without limitation, sample contamination from biological or non-biological sources, specimen marking issues, rare genetic variants interfering with analysis, and other technical issues and limitations.

DISCLAIMERS

This test was developed, and its performance characteristics were determined, by MyOme, Inc., a clinical laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and College of American Pathologist (CAP) accredited to perform high complexity clinical laboratory testing. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA). Test results should always be interpreted by a clinician in the context of clinical and familial data with the availability of genetic counseling when appropriate. MyOme is not responsible for the content of third-party websites referenced in this report.

REVIEWED BY

MyOme

MyOme Example Lab Director

02/13/2025

Date